

Coagulation

# CL 4

**BE**  
Behnk Elektronik  
Diagnostic  
Systems

## The coagulometer for whole blood and blood plasma



**Safely determines all coagulation parameters**

# The optimal measuring station for all coagulation parameters

**The measuring system:** A sensor monitors the ball rotating on the bottom of the cuvette. When a large coagulum is created, the ball is stopped on the wall of the cuvette. A smaller coagulum either slows the ball's movement or directs it towards the centre of the cuvette. The following three parameters lead to the positive detection of a developing coagulum: The ball does not pass the sensor because it has stopped; the ball is no longer detected, because it has been directed toward the centre of the cuvette; the impulse sensor detects a delay in the ball's path. Fibrinolyse measures the time taken for a sample to coagulate before measuring the time it takes for the coagulum to break up.

**The calculator:** The operation is simple and easy to understand. All normal calculation types can be selected in dialogue form with the printer. The programming status can be programmed by using the <E> key.

**The handling:** A register allows up to 40 patient samples to be prepared and then continually processed. Three assigned incubation positions enable a safe and continuous routine. Two reagent mixing positions and an incubation position for a Multipette® round off the instrument. Three positions for reagents and one for the ball dispenser keep the working area organised.

## Features:

Four measuring channels  
Twelve incubation positions  
One incubation position for a Multipette®  
One register for 40 samples  
Three positions for reagent containers  
One position for a ball dispenser  
One quiet thermo-printer  
One RS 232 interface  
Automatically saves calibration curves  
Sample volumes between 150 and 400 µl  
Error control during the coagulation process  
Automatic or manual start of the measurement  
Single or double determination

Operating voltage: 115/230 V ± 10 %  
Frequency range: 50 - 60 Hz  
Power consumption: 60 VA

## Dimensions:

Length: 545 mm  
Width: 300 mm  
Height: 200 mm  
Weight: 13.4 kg

## The basics:

28 Cuvette racks  
250 Balls  
20 Reagent containers Ø 25 mm  
10 Reagent containers Ø 30 mm  
10 Lids for reagent containers Ø 30 mm  
5 Lids for reagent containers Ø 25 mm  
1 Ball dispenser

2 Magnetic stirrers  
1 Multipette® with a start cable  
1 Clean card  
1 Roll of thermo-paper  
10 Combitips® 2.5 ml  
1 Dust cover

Behnk Elektronik GmbH & Co. KG · Hans-Böckler-Ring 27 · D-22851 Norderstedt · Germany

